



Sterilance Lite III

Sterile Safety lancet

Vitrex Sterilance Lite III lancet is intended for blood glucose and other in vitro diagnostic tests. Health workers worldwide face accidental needle-stick injury while caring for patients, putting themselves at risk to bloodborne illnesses including the AIDS virus and Hepatitis.

BUTTON ACTIVATION

USER SAFETY: Needle is fully shielded before and after use.

USER CONVENIENCE: Pre-loaded, pressure activated design.

PATIENT COMFORT: High quality, ultra sharp tri-bevel needle or bi-bevel blade, high speed penetration.

EFFICACY: Variety of point configurations and penetration depths to suit most capillary blood requirements.

PATIENT SAFETY: Single-use only. Sterilised by gamma irradiation.



1



2



3



4



Recommended use

1. Twist off the protective cap.
2. Select fingertip insertion site slightly off centre to minimise any discomfort and maximise blood flow. Hold the lancet firmly against site and press the release button to activate the device.
3. Gently massage the finger in the direction of the puncture site to obtain the required volume of blood.
4. Dispose the safety lancet in a suitable safety container.

WARNING: Do not use, if the protective cap has previously been removed.

Technical Specifications

Colour Code	Specification	Penetration Depth	Article No.	Application
Green ●	18G blade	1.8 mm	373213	Special baby care and maternity units for heel puncture
Orange ●	21G needle	1.8 mm	377213	Haemoglobin (blood banks), cholesterol, blood typing, coagulation, blood gas and other diagnostic tests. For paediatrics and most normal adult skin
Pink ●	21G needle	2.4 mm	378213	As the Orange lancet, but for thicker adult skin
Yellow ●	26G needle	1.8 mm	375213	Diabetes clinics/hospital wards for blood glucose monitoring
Blue ●	26G needle	2.4 mm	376213	As the Yellow lancet, but for thicker adult skin
Purple ●	28G needle	1.8 mm	374213	Diabetes clinics/hospital wards for blood glucose monitoring

Shelf life 5 years

PRN
MEDICAL SUPPLIES

www.prnmed.com.au
sales@prnmed.com.au